

PACO33640

Product Information

Size:

50ug

Reactivity:

Gallus gallus

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May be an autocrine factor that promotes the growth of fibroblasts and is involved in the neoplastic transformation of fibroblasts by v-Src. Chemotactic for peripheral blood mononuclear cells as well as for heterophils.

Gene ID:

CXCL8

Uniprot

P08317

Synonyms:

Interleukin-8 (IL-8) (9E3) (C-X-C motif chemokine 8) (CEF-4) (Chemokine (C-X-C motif ligand 8) (Embryo fibroblast protein 1) (EMF-1), CXCL8, CEF4 EMF1 IL8

Immunogen:

Recombinant Gallus gallus Interleukin-8 protein (17-102AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A