

PACO33696

---

## Product Information

**Size:**

50ug

**Reactivity:**

Salmonella typhi

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Required for the action of colicins K and L and for the stabilization of mating aggregates in conjugation. Serves as a receptor for a number of T-even like phages. Also acts as a porin with low permeability that allows slow penetration of small solutes.

**Gene ID:**

ompA

**Uniprot**

Q8Z7S0

**Synonyms:**

Outer membrane protein A, ompA

**Immunogen:**

Recombinant Salmonella typhi Outer membrane protein A protein (27-349AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A