

PACO33716

Product Information

Size:

50ug

Reactivity:

Escherichia coli

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Facilitates transcription termination by a mechanism that involves rho binding to the nascent RNA, activation of rho's RNA-dependent ATPase activity, and release of the mRNA from the DNA template. RNA-dependent NTPase which utilizes all four ribonucleoside triphosphates as substrates.

Gene ID:

rho

Uniprot

P0AG30

Synonyms:

Transcription termination factor Rho (EC 3.6.4. -) (ATP-dependent helicase Rho), rho, nitA psuA rnsC sbaA tsu

Immunogen:

Recombinant Escherichia coli Transcription termination factor Rho protein (1-419AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A