

# IL18BP Antibody



PACO33740

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

Isoform A binds to IL-18 and inhibits its activity. Functions as an inhibitor of the early TH1 cytokine response.

**Gene ID:**

IL18BP

**Uniprot**

O95998

**Synonyms:**

Interleukin-18-binding protein (IL-18BP) (Tadekinig-alfa), IL18BP

**Immunogen:**


Recombinant Human Interleukin-18-binding protein (31-183AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human placenta tissue using PACO33740 at dilution of 1:100.