

IL18BP Antibody, HRP conjugated



PACO33741

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Isoform A binds to IL-18 and inhibits its activity. Functions as an inhibitor of the early TH1 cytokine response.

Gene ID:

IL18BP

Uniprot

O95998

Synonyms:

Interleukin-18-binding protein (IL-18BP) (Tadekinig-alfa), IL18BP

Immunogen:

Recombinant Human Interleukin-18-binding protein (31-183AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A