

# Il18bp Antibody



PACO33744

---

## Product Information

**Size:**

50ug

**Reactivity:**

Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:10-1:100

**Protein Background:**

Binds to IL-18 and inhibits its activity. Functions as an inhibitor of the early TH1 cytokine response.

**Gene ID:**

Il18bp

**Uniprot**

Q9Z0M9

**Synonyms:**

Interleukin-18-binding protein (IL-18BP) (Interferon gamma-inducing factor-binding protein), Il18bp, Igifbp

**Immunogen:**


Recombinant Mouse Interleukin-18-binding protein (29-193AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human colon cancer using PACO33744 at dilution of 1:20.