

Il18bp Antibody, HRP conjugated



PACO33745

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Binds to IL-18 and inhibits its activity. Functions as an inhibitor of the early TH1 cytokine response.

Gene ID:

Il18bp

Uniprot

Q9Z0M9

Synonyms:

Interleukin-18-binding protein (IL-18BP) (Interferon gamma-inducing factor-binding protein), Il18bp, Igifbp

Immunogen:

Recombinant Mouse Interleukin-18-binding protein (29-193AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A