

PACO33760

---

## Product Information

**Size:**

50ug

**Reactivity:**

Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

This isozyme may play a role in skeletal mineralization.

**Gene ID:**

Alpl

**Uniprot**

P09242

**Synonyms:**

Alkaline phosphatase, tissue-nonspecific isozyme (AP-TNAP) (TNSALP) (EC 3.1.3.1) (Alkaline phosphatase 2) (Alkaline phosphatase liver/bone/kidney isozyme), Alpl, Akp-2 Akp2

**Immunogen:**

Recombinant Mouse Alkaline phosphatase, tissue-nonspecific isozyme protein (20-503AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A