

Alpl Antibody, Biotin conjugated



PACO33763

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

This isozyme may play a role in skeletal mineralization.

Gene ID:

Alpl

Uniprot

P09242

Synonyms:

Alkaline phosphatase, tissue-nonspecific isozyme (AP-TNAP) (TNSALP) (EC 3.1.3.1) (Alkaline phosphatase 2) (Alkaline phosphatase liver/bone/kidney isozyme), Alpl, Akp-2 Akp2

Immunogen:

Recombinant Mouse Alkaline phosphatase, tissue-nonspecific isozyme protein (20-503AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A