

Muc5ac Antibody



PACO33780

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

Muc5ac

Uniprot

E9PWB6

Synonyms:

Mucin 5, subtypes A and C, tracheobronchial/gastric, Muc5ac

Immunogen:

Recombinant Mouse Mucin 5, subtypes A and C, tracheobronchial/gastric protein (2452-2721AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A