

MT1F Antibody



PACO33816

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Metallothioneins have a high content of cysteine residues that bind various heavy metals; these proteins are transcriptionally regulated by both heavy metals and glucocorticoids.

Gene ID:

MT1F

Uniprot

P04733

Synonyms:

Metallothionein-1F (MT-1F) (Metallothionein-IF) (MT-IF), MT1F

Immunogen:

Recombinant Human Metallothionein-1F protein (1-61AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A