

# MT1M Antibody



PACO33824

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Metallothioneins have a high content of cysteine residues that bind various heavy metals; these proteins are transcriptionally regulated by both heavy metals and glucocorticoids.

**Gene ID:**

MT1M

**Uniprot**

Q8N339

**Synonyms:**

Metallothionein-1M (MT-1M) (Metallothionein-IM) (MT-IM), MT1M, MT1K

**Immunogen:**

Recombinant Human Metallothionein-1M protein (1-61AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A