

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Mediates zinc uptake. May function as a major endogenous zinc uptake transporter in many cells of the body. Responsible for the rapid uptake and accumulation of physiologically effective zinc in prostate cells.

Gene ID:

SLC39A1

Uniprot

Q9NY26

Synonyms:

Zinc transporter ZIP1 (Solute carrier family 39 member 1) (Zinc-iron-regulated transporter-like) (Zrt- and Irt-like protein 1) (ZIP-1) (hZIP1), SLC39A1, IRT1 ZIP1 ZIRTL

Immunogen:

Recombinant Human Zinc transporter ZIP1 protein (126-179AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A