

# poxB Antibody



PACO33880

---

## Product Information

**Size:**

50ug

**Reactivity:**

Escherichia coli

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

poxB

**Uniprot**

P07003

**Synonyms:**

Pyruvate dehydrogenase [ubiquinone] (EC 1.2.5.1) (Pyruvate oxidase) (POX) [Cleaved into: Alpha-peptide], poxB

**Immunogen:**

Recombinant Escherichia coli Pyruvate dehydrogenase [ubiquinone] protein (1-572AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A