

PACO33896

---

## Product Information

**Size:**

50ug

**Reactivity:**

Pig

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Pulmonary surfactant-associated proteins promote alveolar stability by lowering the surface tension at the air-Liquid interface in the peripheral air spaces. SP-B increases the collapse pressure of palmitic acid, to nearly 70 millinewtons per meter.

**Gene ID:**

SFTPB

**Uniprot**

P15782

**Synonyms:**

Pulmonary surfactant-associated protein B (SP-B) (8 kDa protein) (Pulmonary surfactant-associated proteolipid SPL(Phe)), SFTPB, SFTP3

**Immunogen:**

Recombinant Pig Pulmonary surfactant-associated protein B protein (1-79AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A