

COL17A1 Antibody, FITC conjugated



PACO33902

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May play a role in the integrity of hemidesmosome and the attachment of basal keratinocytes to the underlying basement membrane.

Gene ID:

COL17A1

Uniprot

Q9UMD9

Synonyms:

Collagen alpha-1(XVII) chain (180 kDa bullous pemphigoid antigen 2) (Bullous pemphigoid antigen 2) [Cleaved into: 120 kDa linear IgA disease antigen (120 kDa linear IgA dermatosis antigen) (Linear IgA disease antigen 1) (LAD-1); 97 kDa linear IgA disease antigen (97 kDa linear IgA bullous dermatosis antigen) (97 kDa LAD antigen) (97-LAD) (Linear IgA bullous disease antigen of 97 kDa) (LABD97)], COL17A1, BP180 BPAG2

Immunogen:

Recombinant Human Collagen alpha-1(XVII) chain protein (1253-1497AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A