

Tubulin beta-1 chain Antibody



PACO33928

Product Information

Size:

50ug

Reactivity:

Gallus gallus

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Tubulin is the major constituent of microtubules. It binds two moles of GTP, one at an exchangeable site on the beta chain and one at a non-exchangeable site on the alpha-chain.

Gene ID:

Uniprot

P09203

Synonyms:

Tubulin beta-1 chain (Beta-tubulin class-I)

Immunogen:

Recombinant Gallus gallus Tubulin beta-1 chain protein (1-439AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A