

Serine protease inhibitor Kazal-type 2-like protein Antibody
PACO33968



Product Information

Size:

50ug

Reactivity:

Cyprinus carpio

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

Uniprot

A0A059T4U4

Synonyms:

Immunogen:

Recombinant Cyprinus carpio Serine protease inhibitor Kazal-type 2-like protein (18-77AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A