

HLA-DRB1 Antibody, FITC conjugated



PACO33974

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

HLA-DRB1

Uniprot

Q95IE3

Synonyms:

HLA class II histocompatibility antigen, DRB1-12 beta chain (MHC class II antigen DRB1*12) (DR-12) (DR12), HLA-DRB1

Immunogen:

Recombinant Human HLA class II histocompatibility antigen, DRB1-12 & beta; chain protein (30-266AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A