

PACO33984

Product Information

Size:

50ug

Reactivity:

Enterobacteria phage H19B

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

The B subunit is responsible for the binding of the holotoxin to specific receptors on the target cell surface, such as globotriaosylceramide (Gb3) in human intestinal microvilli.

Gene ID:

stxB

Uniprot

P69179

Synonyms:

Shiga-like toxin 1 subunit B (SLT-1 B subunit) (SLT-1b) (SLT-1b) (Verocytotoxin 1 subunit B) (Verotoxin 1 subunit B), stxB, sltB

Immunogen:

Recombinant Enterobacteria phage H19B Shiga-like toxin 1 subunit B protein (21-89AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A