

glsA1 Antibody



PACO34000

Product Information

Size:

50ug

Reactivity:

Escherichia coli

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000

Protein Background:

Gene ID:

glsA1

Uniprot

P77454

Synonyms:

Glutaminase 1 (EC 3.5.1.2), glsA1, ybaS


Immunogen:

Recombinant Escherichia coli Glutaminase 1 protein (1-310AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot. All lanes: glsA1 antibody at 2 μ g/ml + DH5a whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 33 kDa. Observed band size: 33 kDa.