

Pectate lyase 1 Antibody



PACO34032

Product Information

Size:

50ug

Reactivity:

Juniperus ashei

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Has low pectate lyase activity.

Gene ID:

Uniprot

P81294

Synonyms:

Pectate lyase 1 (EC 4.2.2.2) (Major pollen allergen Jun a 1) (allergen Jun a 1)

Immunogen:

Recombinant Juniperus ashei Pectate lyase 1 protein (22-367AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A