

**Ribosome-inactivating protein PD-L1/PD-L2
Antibody, Biotin conjugated
PACO34039**



Product Information

Size:

50ug

Reactivity:

Phytolacca dioica

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Inhibits protein synthesis. Cleaves supercoiled pBR322 dsDNA.

Gene ID:

Uniprot

P84853

Synonyms:

Ribosome-inactivating protein PD-L1/PD-L2 (EC 3.2.2.22) (rRNA N-glycosidase PD-L1/PD-L2)

Immunogen:

Recombinant Phytolacca dioica Ribosome-inactivating protein PD-L1/PD-L2 protein (1-261AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A