

Alpha-amanitin proprotein Antibody



PACO34040

Product Information

Size:

50ug

Reactivity:

Amanita phalloides

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Major toxin from Amanita phalloideae. Acts by binding non-competitively to RNA polymerase II and greatly slowing the elongation of transcripts from target promoters.

Gene ID:**Uniprot**

P85421

Synonyms:

Alpha-amanitin proprotein (Alpha-amanitin) (Amatoxin) (Gamma-amanitin)

Immunogen:

Peptide sequence from Amanita phalloides Alpha-amanitin proprotein protein.

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A