

**Ribosome-inactivating protein lychnin Antibody,
HRP conjugated**
PACO34047



Product Information

Size:

50ug

Reactivity:

Silene chaledonica

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Ribosome-inactivating protein of type 1, inhibits protein synthesis in animal cells. Inhibits cell-free translation in rabbit reticulocyte lysate system with an IC(50) of 0.17 nM.

Gene ID:

Uniprot

P85101

Synonyms:

Ribosome-inactivating protein lychnin (EC 3.2.2.22)

Immunogen:

Recombinant Silene chaledonica Ribosome-inactivating protein lychnin protein (1-234AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A