

Major pollen allergen Dac g 4 Antibody



PACO34050

Product Information

Size:

50ug

Reactivity:

Dactylis glomeRata

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

Uniprot

P82946

Synonyms:

Major pollen allergen Dac g 4 (allergen Dac g 4) (Fragments)

Immunogen:

Recombinant Dactylis glomerata Major pollen allergen Dac g 4 protein (1-55AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A