

Major pollen allergen Dac g 4 Antibody, FITC conjugated
PACO34052



Product Information

Size:

50ug

Protein Background:

Reactivity:

Dactylis glomerata

Gene ID:

Source:

Rabbit

Uniprot

P82946

Isotype:

IgG

Synonyms:

Major pollen allergen Dac g 4 (allergen Dac g 4) (Fragments)

Applications:

ELISA

Immunogen:

Recombinant Dactylis glomerata Major pollen allergen Dac g 4 protein (1-55AA).

Recommended dilutions:

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A