

Ribosome-inactivating protein charybdin
Antibody, Biotin conjugated
PACO34057



Product Information

Size:

50ug

Reactivity:

Drimia maritima

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Inhibits translation in rabbit reticulocytes.

Gene ID:

Uniprot

P84786

Synonyms:

Ribosome-inactivating protein charybdin (EC 3.2.2.22) (rRNA N-glycosidase)

Immunogen:

Recombinant Drimia maritima Ribosome-inactivating protein charybdin protein (1-257AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A