

Pollen allergen Dac g 3 Antibody



PACO34062

Product Information

Size:

50ug

Reactivity:

Dactylis glomeRata

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

Uniprot

P93124

Synonyms:

Pollen allergen Dac g 3 (Allergen Dac g III) (allergen Dac g 3) (Fragment)

Immunogen:

Recombinant Dactylis glomerata Pollen allergen Dac g 3 protein (1-96AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A