

**Pollen allergen Dac g 3 Antibody, FITC conjugated**  
**PACO34064**



---

**Product Information**

**Size:**

50ug

**Protein Background:**

**Reactivity:**

Dactylis glomeRata

**Gene ID:**

**Source:**

Rabbit

**Uniprot**

P93124

**Isotype:**

IgG

**Synonyms:**

Pollen allergen Dac g 3 (Allergen Dac g III) (allergen Dac g 3) (Fragment)

**Applications:**

ELISA

**Immunogen:**

Recombinant Dactylis glomerata Pollen allergen Dac g 3 protein (1-96AA).

**Recommended dilutions:**

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A