

Major allergen Api g 1, isoallergen 2 Antibody



PACO34070

Product Information

Size:

50ug

Reactivity:

Apium graveolens

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

Uniprot

P92918

Synonyms:

Major allergen Api g 1, isoallergen 2 (Allergen Api g 1.0201) (allergen Api g 1)

Immunogen:

Recombinant Apium graveolens Major allergen Api g 1, isoallergen 2 protein (1-159AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A