

# F Antibody, FITC conjugated



PACO34136

---

## Product Information

**Size:**

50ug

**Reactivity:**

Canine distemper virus

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

F

**Uniprot**

P12569

**Synonyms:**

Fusion glycoprotein F0 [Cleaved into: Fusion glycoprotein F2; Fusion glycoprotein F1], F

**Immunogen:**

Recombinant Canine distemper virus Fusion glycoprotein F0 protein (136-608AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A