

MDR1 Antibody



PACO34138

Product Information

Size:

50ug

Reactivity:

Plasmodium falciparum

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Energy-dependent efflux pump responsible for decreased drug accumulation in multidrug-resistant cells.

Gene ID:

MDR1

Uniprot

P13568

Synonyms:

Multidrug resistance protein (EC 3.6.3.44) (Chloroquine resistance protein), MDR1

Immunogen:

Recombinant Plasmodium falciparum Multidrug resistance protein (1084-1419AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A