

# fla Antibody



PACO34170

---

## Product Information

**Size:**

50ug

**Reactivity:**

Borrelia burgdorferi

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Component of the core of the flagella.

**Gene ID:**

fla

**Uniprot**

P11089

**Synonyms:**

Flagellar filament 41 kDa core protein (Flagellin) (41 kDa antigen) (P41), fla

**Immunogen:**

Recombinant Borrelia burgdorferi Flagellar filament 41 kDa core protein (1-336AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A