

PYC1 Antibody



PACO34174

Product Information

Size:

50ug

Reactivity:

Baker's yeast

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Pyruvate carboxylase catalyzes a 2-step reaction, involving the ATP-dependent carboxylation of the covalently attached biotin in the first step and the transfer of the carboxyl group to pyruvate in the second.

Gene ID:

PYC1

Uniprot

P11154

Synonyms:

Pyruvate carboxylase 1 (EC 6.4.1.1) (Pyruvic carboxylase 1) (PCB 1), PYC1, PYV

Immunogen:

Recombinant Baker's yeast Pyruvate carboxylase 1 protein (18-470AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A