

Outer capsid protein VP4 Antibody, Biotin conjugated

PACO34201



Product Information

Size:

50ug

Reactivity:

Rotavirus A

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:****Uniprot**

P11200

Synonyms:

Outer capsid protein VP4 (Hemagglutinin) [Cleaved into: Outer capsid protein VP8*; Outer capsid protein VP5*]

Immunogen:

Recombinant Rotavirus A Outer capsid protein VP4 protein (247-479AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A