

Pla2g2a Antibody



PACO34222

Product Information

Size:

50ug

Reactivity:

Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Thought to participate in the regulation of the phospholipid metabolism in biomembranes including eicosanoid biosynthesis. Catalyzes the calcium-dependent hydrolysis of the 2-acyl groups in 3-sn-phosphoglycerides.

Gene ID:

Pla2g2a

Uniprot

P14423

Synonyms:

Phospholipase A2, membrane associated (EC 3.1.1.4) (GLIC sPLA2) (Group IIA phospholipase A2) (Phosphatidylcholine 2-acylhydrolase 2A), Pla2g2a

Immunogen:

Recombinant Rat Phospholipase A2, membrane associated protein (22-146AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A