

toxB Antibody



PACO34238

Product Information

Size:

50ug

Reactivity:

Peptoclostridium difficile

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

toxB

Uniprot

P18177

Synonyms:

Toxin B (EC 3.4.22. -), toxB, tcdB

Immunogen:

Recombinant Peptoclostridium difficile Toxin B protein (1751-2366AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A