

# Sialic acid-binding lectin Antibody



PACO34246

---

## Product Information

**Size:**

50ug

**Reactivity:**

Rana japonica

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

The S-lectins in frog eggs may be involved in the fertilization and development of the frog embryo. This lectin preferentially agglutinates a large variety of tumor cells, but it does not agglutinate non-transformed cells and erythrocytes.

**Gene ID:**

**Uniprot**

P18839

**Synonyms:**

Sialic acid, binding lectin (EC 3.1.27. -)

**Immunogen:**

Recombinant Rana japonica Sialic acid, binding lectin protein (1-111AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A