

# Major pepsin inhibitor 3 Antibody, HRP conjugated

PACO34255



---

## Product Information

**Size:**

50ug

**Reactivity:**

Ascaris suum

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

This is an inhibitor of the aspartic protease pepsin.

**Gene ID:****Uniprot**

P19400

**Synonyms:**

Major pepsin inhibitor 3 (PI-3)

**Immunogen:**

Recombinant Ascaris suum Major pepsin inhibitor 3 protein (21-169AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A