

Major pepsin inhibitor 3 Antibody, Biotin conjugated

PACO34257



Product Information

Size:

50ug

Reactivity:

Ascaris suum

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

This is an inhibitor of the aspartic protease pepsin.

Gene ID:**Uniprot**

P19400

Synonyms:

Major pepsin inhibitor 3 (PI-3)

Immunogen:

Recombinant Ascaris suum Major pepsin inhibitor 3 protein (21-169AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A