

IFNT1 Antibody



PACO34262

Product Information

Size:

50ug

Reactivity:

Bovine

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

IFNT1

Uniprot

P15696

Synonyms:

Interferon tau-1 (IFN-tau-1) (Antiluteolysin) (Trophoblast antiluteolytic protein) (Trophoblast protein 1) (TP-1) (Trophoblastin), IFNT1

Immunogen:

Recombinant Bovine Interferon tau-1 protein (24-195AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A