

# P Antibody



PACO34294

---

## Product Information

**Size:**

50ug

**Reactivity:**

Hepatitis B virus genotype A2 subtype adw2

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:****Gene ID:**

P

**Uniprot**

P17100

**Synonyms:**

Protein P [Includes: DNA-directed DNA polymerase (EC 2.7.7.7); RNA-directed DNA polymerase (EC 2.7.7.49); Ribonuclease H (EC 3.1.26.4)], P

**Immunogen:**

Recombinant Hepatitis B virus genotype A2 subtype adw2 Protein P protein (349-692AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A