

Klk1b22 Antibody



PACO34314

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Glandular kallikreins cleave Met-Lys and Arg-Ser bonds in kininogen to release Lys-bradykinin.

Gene ID:

Klk1b22

Uniprot

P15948

Synonyms:

Kallikrein 1-related peptidase b22 (EC 3.4.21.35) (Beta-NGF-endorpeptidase) (Epidermal growth factor-binding protein type A) (EGF-BP A) (Glandular kallikrein K22) (mGK-22) (Nerve growth factor beta chain endopeptidase) (Tissue kallikrein 22), Klk1b22, Klk-22, Klk22

Immunogen:

Recombinant Mouse Kallikrein 1-related peptidase b22 protein (25-259AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A