

# KP6 killer toxin Antibody



PACO34326

---

## Product Information

**Size:**

50ug

**Reactivity:**

Ustilago maydis P6 virus

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

This protein is lethal to sensitive cells of the same or related species. The KP6 alpha subunit is known to recognize some cellular receptors before interaction of the complex with KP6 beta, precipitating cell death.

**Gene ID:****Uniprot**

P16948

**Synonyms:**

KP6 killer toxin (Killer protein 6) [Cleaved into: KP6 killer toxin subunit alpha (VP10); KP6 killer toxin subunit beta (VP12.5)]

**Immunogen:**

Recombinant Ustilago maydis P6 virus KP6 killer toxin protein (28-105AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A