

KP6 killer toxin Antibody, HRP conjugated



PACO34327

Product Information

Size:

50ug

Reactivity:

Ustilago maydis P6 virus

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

This protein is lethal to sensitive cells of the same or related species. The KP6 alpha subunit is known to recognize some cellular receptors before interaction of the complex with KP6 beta, precipitating cell death.

Gene ID:**Uniprot**

P16948

Synonyms:

KP6 killer toxin (Killer protein 6) [Cleaved into: KP6 killer toxin subunit alpha (VP10); KP6 killer toxin subunit beta (VP12.5)]

Immunogen:

Recombinant Ustilago maydis P6 virus KP6 killer toxin protein (28-105AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A