

# Ribosome-inactivating protein karasurin-C Antibody PACO34330



---

## Product Information

**Size:**

50ug

**Reactivity:**

Trichosanthes kirilowii

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Abortion-inducing protein. It inactivates eukaryotic 60S ribosomal subunits.

**Gene ID:****Uniprot**

P24478

**Synonyms:**

Ribosome-inactivating protein karasurin-C (EC 3.2.2.22) (rRNA N-glycosidase) [Cleaved into: Ribosome-inactivating protein karasurin-A]

**Immunogen:**

Recombinant Trichosanthes kirilowii Ribosome-inactivating protein karasurin-C protein (22-270AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A