

**Ribosome-inactivating protein karasurin-C
Antibody, FITC conjugated**
PACO34332



Product Information

Size:

50ug

Reactivity:

Trichosanthes kirilowii

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Abortion-inducing protein. It inactivates eukaryotic 60S ribosomal subunits.

Gene ID:

Uniprot

P24478

Synonyms:

Ribosome-inactivating protein karasurin-C (EC 3.2.2.22) (rRNA N-glycosidase) [Cleaved into: Ribosome-inactivating protein karasurin-A]

Immunogen:

Recombinant Trichosanthes kirilowii Ribosome-inactivating protein karasurin-C protein (22-270AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A