

ampC Antibody



PACO34346

Product Information

Size:

50ug

Reactivity:

Pseudomonas aeruginosa

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

ampC

Uniprot

P24735

Synonyms:

Beta-lactamase (EC 3.5.2.6) (Cephalosporinase), ampC

Immunogen:

Recombinant Pseudomonas aeruginosa Beta-lactamase protein (27-397AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A