

Reg3b Antibody, HRP conjugated



PACO34351

Product Information

Size:

50ug

Reactivity:

Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Bactericidal C-type lectin which acts against several intestinal Gram-positive bacteria and Gram-negative bacteria. Lacks antibacterial activity against *S. typhimurium*. May play a role in protection against infection with *S. enteritidis* by inhibiting its translocation from the gut lumen into intestinal tissues and further extraintestinal tissues.

Gene ID:

Reg3b

Uniprot

P25031

Synonyms:

Regenerating islet-derived protein 3-beta (REG-3-beta) (Pancreatitis-associated protein 1) (Peptide 23) (REG-2) (Regenerating islet-derived protein III-beta) (Reg III-beta) [Cleaved into: Regenerating islet-derived protein 3-beta 16.5 kDa form; Regenerating islet-derived protein 3-beta 15 kDa form], Reg3b, Pap Pap1 Reg2

Immunogen:

Recombinant Rat Regenerating islet-derived protein 3-β protein (27-175AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A